



Guide 5

Personal Voice and voice banking

Keep your own voice, and set exactly how Rejoin speaks.

If speech is still clear, recording a voice now is one of the most valuable things to do. Rejoin can then speak in that voice for years to come. This guide covers banking a voice and tuning how Rejoin sounds.


For patients and caregivers. Worth doing early.

1. Do it early

For a progressive condition, bank a voice while speech is still strong. It takes about fifteen minutes and only has to be done once. Even if speech is already changing, it is still worth doing.

2. Record a Personal Voice


Personal Voice is an Apple feature, recorded in the device's Settings: Accessibility, then Personal Voice, then Create a Personal Voice. You read a set of phrases aloud and the device builds the voice, entirely on the device.



Voices in Rejoin.

3. Use it in Rejoin

Once it is recorded, open Rejoin, go to Me, then Voice and speech, and choose your Personal Voice. From then on Rejoin speaks in it.



Voice and speech.

4. A cloned voice with Rejoin+

With Rejoin+, Rejoin can also create a natural cloned voice from a short recording, a good option if a Personal Voice was never banked. Either way, the goal is the same: your words in a voice that sounds like you.

5. Tune how it speaks

In Voice and speech you can set the speaking rate and pitch, hear a sample, and teach Rejoin how to say tricky names and words under Say it your way, so a grandchild's name or a hometown always comes out right.